


Product Name: Recombinant Human MMP14 Protein
Catalog #: PEH4485



Summary

Name	MMP14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P50281
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LPDDRRGIQ QLYGGESGFP TKMPPQRRT SRPSVPDKPK NPTYGPNICD GNFDTVAMLR GEMFVFKERW FWRVRNNQVM DGYPMPIGQF WRGLPASINT AYERKDGKQV FFKGDKHWVF DEASLEPGYP KHIKELGRGL PTDKIDAALF WMPNGKTYFF RGNKYRFNE ELRAVDSEYP KNIKVWEGIP ESPRGSFMGS DEVFTYFYKG NKYWKFNQK LKVEPGYPKS ALRDWMGCPS GGRPDEGTEE ETEVIIEVD EEGGGAVSAA

image



Background

Alternative Names	Matrix metalloproteinase-14 (MMP-14) (EC 3.4.24.80) (MMP-X1) (Membrane-type matrix metalloproteinase 1) (MT-MMP 1) (MTMMP1) (Membrane-type-1 matrix metalloproteinase) (MT1-MMP) (MT1MMP)
--------------------------	---

Product Name: Recombinant Human MMP14 Protein
Catalog #: PEH4485



Note

For Research Use Only , Not for Diagnostic Use