



**Product Name: Recombinant Human MMP14 Protein**  
**Catalog #: PEH4486**



## Summary

<b>Name</b>	MMP14_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-SUMO & C-His
<b>Accession #</b>	P50281
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	44 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LPDDRRGIQ QLYGGESGFP TKMPPQRRT SRPSVPDKPK NPTYGPNICD GNFDTVAMLR GEMFVFKERW FWRVRNNQVM DGYPMPIGQF WRGLPASINT AYERKDGKQV FFKGDKHWVF DEASLEPGYP KHIKELGRGL PTDKIDAALF WMPNGKTYFF RGNKYRFNE ELRAVDSEYP KNIKVWEGIP ESPRGSFMGS DEVFTYFYKG NKYWKFNQK LKVEPGYPS ALRDWMGCPS GGRPDEGTEE ETEVIIEVD EEGGAVSAA

## image



## Background

**Alternative Names** Matrix metalloproteinase-14 (MMP-14) (EC 3.4.24.80) (MMP-X1) (Membrane-type matrix metalloproteinase 1) (MT-MMP 1) (MTMMP1) (Membrane-type-1 matrix metalloproteinase) (MT1-MMP) (MT1MMP)

**Product Name: Recombinant Human MMP14 Protein**  
**Catalog #: PEH4486**



---

**Note**

For Research Use Only , Not for Diagnostic Use