


Product Name: Recombinant Human DUSP3 Protein
Catalog #: PEH4488



Summary

Name	DUS3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P51452
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSGSFELSVQ DLNDLLSDGS GCYSLPSQPC NEVTPRIYVG NASVAQDIPK LQKLGITHVL NAAEGRSFMH VNTNANFYKD SGITYLGIKA NDTQEFNLSA YFERAADFID QALAQKNGRV LVHCREGYSR SPTLVIAYLM MRQKMDVKSA LSIVRQNREI GPNDGFLAQL CQLNDRLAKE GKLKP

image



Background

Alternative Names	Dual specificity protein phosphatase 3 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase VHR) (Vaccinia H1-related phosphatase) (VHR)
--------------------------	---

Note

Product Name: Recombinant Human DUSP3 Protein
Catalog #: PEH4488



For Research Use Only , Not for Diagnostic Use