


Product Name: Recombinant Human PRELP Protein
Catalog #: PEH4489



Summary

Name	PRELP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P51888
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPGRRPRPRP RPTPSFPQPD EPAEPTDLPP PLPPGPPSIF PDCPRECYCP PDFPSALYCD SRNLRKVPVI PPRIHYLYLQ NNFITELPVE SFQNATGLRW INLDNNRIRK IDQRVLEKLP GLVFLYMEKN QLEEVPSALP RNLEQLRLSQ NHISRIPPGV FSKLENLLLL DLQHNRLSDG VFKPDTFHGL KNLMQLNLAH NILRKMPPRV PTAIHQLYLD SNKIETIPNG YFKSFPNLAF IRLNYNKLTD

image



Background

Alternative Names Prolargin (Proline-arginine-rich end leucine-rich repeat protein)

Note

Product Name: Recombinant Human PRELP Protein
Catalog #: PEH4489



For Research Use Only , Not for Diagnostic Use