



Product Name: Recombinant Human EFNB2 Protein
Catalog #: PEH4492



Summary

Name	EFNB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P52799
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EPIYWSSNS KFLPGQGLVL YPQIGDKLDI ICPKVDSKTV GQYEYKQVYM VDKDQADRCT IKKENTPLLN CAKPDQDIKF TIKFQEFSPN LWGLEFQKNK DYYIISTSG SLEGLDNQEG GVCQTRAMKI

image



Background

Alternative Names	Ephrin-B2 (EPH-related receptor tyrosine kinase ligand 5) (LERK-5) (HTK ligand) (HTK-L)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use