


Product Name: Recombinant Human METAP1 Protein
Catalog #: PEH4493



Summary

Name	MAP11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P53582
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KREVSSWTVE GDINTDPWAG YRYTGKLRPH YPLMPTRPVP SYIQRPDYAD HPLGMSESEQ ALKGTSQLKL LSSEDIEGMR LVCRLAREVL DVAAGMIKPG VTTEEIDHAV HLACIARNCY PSPLNYYNFP KSCCTSVNEV ICHGIPDRRP LQEGDIVNVD ITLYRNGYHG DLNETFFVGE VDDGARKLVQ TTYECLMQAI DAVKPGVRYR ELGNIIQKHA QANGFSVVRS YCGHGIHKL F HTAPNVPHYA KNKAVGVMKS GHVFTIEPMI

image



Background

Alternative Names	Methionine aminopeptidase 1 (MAP 1) (MetAP 1) (EC 3.4.11.18) (Peptidase M 1)
--------------------------	--

Product Name: Recombinant Human METAP1 Protein
Catalog #: PEH4493



Note

For Research Use Only , Not for Diagnostic Use