


Product Name: Recombinant Human EPHB4 Protein
Catalog #: PEH4497



Summary

Name	EPHB4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P54760
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ARTDPRGAPC TTPPSAPRSV VSRLNGSSLH LEWSAPLESG GREDLTYALR CRECRPGGSC APCGGDLTFD PGPRDLVEPW VVVRGLRPDF TYTFEVTALN GVSSLATGPV PFEPVNVTTD REVPPAVSDI RVTRSSPSSL SLAWAVPRAP SGAVLDYEVK YHEKGAEGPS SVRFLKTSN RAELRGLKRG ASYLVQVRAR SEAGYGPFQ EHHSTQLDE

image



Background

Alternative Names Ephrin type-B receptor 4 (EC 2.7.10.1) (Hepatoma transmembrane kinase) (Tyrosine-protein kinase TYRO11)

Note

Product Name: Recombinant Human EPHB4 Protein
Catalog #: PEH4497



For Research Use Only , Not for Diagnostic Use