


Product Name: Recombinant Human PLTP Protein
Catalog #: PEH4498



Summary

Name	PLTP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P55058
Host	E.coli
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GCKIRVTSKA LELVKQEGLR FLEQELETIT IPDLRGKEGH FYYNISEVKV TELQLTSSSEL DFQPQQELML QITNASLGLR FRRQLLYWFF YDGGYINASA EGVSIRTGLE LSRDPAGRMK VSNVSCQASV SRMHAAFGGT FKKVYDFLST FITSGMRFLN NQQICPVLYH AGTVLLNSLL DTVPVRSSVD ELVGIDYSLM KDPVASTSNL DMDFRGAFFP LTERNWSLPN RAVEPQLQEE ERMVYVAFSE FFDSDAMESY FRAGALQLLL VGDKVPHDLD MLLRATYFGS IVLLSPAVID SPLKLELRVL APPRCTIKPS GTTISVTASV TIALVPPDQP

image



Background

Alternative Names Phospholipid transfer protein (Lipid transfer protein II)

Product Name: Recombinant Human PLTP Protein
Catalog #: PEH4498



Note

For Research Use Only , Not for Diagnostic Use