


Product Name: Recombinant Human CASP9 Protein
Catalog #: PEH4499



Summary

Name	CASP9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His
Accession #	P55211
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDEADRRLLR RCRLRLVEEL QVDQLWDALL SRELFPHMI EDIQRAGSGS RRDQARQLII DLETRGSQAL PLFISCLEDT GQDMLASFLR TNRQAAKLSK PTLENLTPVV LRPEIRKPEV LRPETPRPVD IGSGGFQDVG ALES LRGNAD LAYILSMEPC GHCLINNPNV FCRESGLRTR TGSNIDCEKL RRRFSSLHFM VEVKGDLTAK KMLVLALELA QQDHGALDCC VVVILSHGCQ ASHLQFPGAV YGTGDCPVSV EKIVNIFNGT SCPSLGKPK LFFIQACGGE QKDHGFEVAS TSPEDSPGS NPEPD

image



Background

Alternative Names	Caspase-9 (CASP-9) (EC 3.4.22.62) (Apoptotic protease Mch-6) (Apoptotic protease-activating factor 3) (APAF-3) (ICE-like apoptotic protease 6) (ICE-
--------------------------	--

Product Name: Recombinant Human CASP9 Protein
Catalog #: PEH4499



LAP6) [Cleaved into: Caspase-9 subunit p35; Caspase-9 subunit p10]

Note

For Research Use Only , Not for Diagnostic Use