



Product Name: Recombinant Human FOXA1 Protein
Catalog #: PEH4501



Summary

Name	FOXA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P55317
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QQRWQNSIRH SLSFNDCFVK VARSPDKPGK GSYWTLHPDS GNMFENGCYL RRQKRFKCEK QPGAGGGGGS GSGGSGAKGG PESRKDPSGA SNPSADSPLH RGVHGKTGQL EGAPAPGPAA SPQTLDHSGA TATGGASELK TPASSTAPPI SSGPGALASV PASHPAHGLA PHESQLHLKG DPHYSFNHPF SINNLMSSE QQHKLDFKAY EQALQYSPYG

image



Background

Alternative Names Hepatocyte nuclear factor 3-alpha (HNF-3-alpha) (HNF-3A) (Forkhead box protein A1) (Transcription factor 3A) (TCF-3A)

Note

Product Name: Recombinant Human FOXA1 Protein
Catalog #: PEH4501



For Research Use Only , Not for Diagnostic Use