


Product Name: Recombinant Human FAM3B Protein
Catalog #: PEH4504



Summary

Name	FAM3B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & N-His
Accession #	P58499
Host	E.coli
Species	Human
Molecular Mass	55 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGDNSGPMTK FIQSAAPKSL LFMVTYDDGS TRLNNDAKNA IEALGSKEIR NMKFRSSWVF IAAKGLELPS EIQRKINHS DAKNNRYSGW PAEIQIEGCI

image



Background

Alternative Names	Protein FAM3B (Cytokine-like protein 2-21) (Pancreatic-derived factor) (PANDER)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use