


Product Name: Recombinant Human DLG4 Protein
Catalog #: PEH4506



Summary

Name	DLG4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P78352
Host	E.coli
Species	Human
Molecular Mass	64 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IKLIKGPKGL GFSIAGGVGN QHIPGNSIY VTKIIEGGAA HKDGRQLIGD KILAVNSVGL EDVMHEDAVA ALKNTYDVVY LKVAKPSNAY LSDSYAPPDI TTSYSQHLDN EISHSSYLGT DYPTAMTPTS PRRYSPVAKD LLGEEDIPRE PRRIVIHRS TGLGFNIVGG EDGEGIFISF ILAGGPADLS GELRKGQIL SVNGVDLRNA SHEQAAIALK NAGQVTIIA QYKPEEYSRF EAKIHDREQ LMNSSLGSGT ASLRSNPKRG FYIRALFDYD KTKDCGFLSQ ALSFRFGDVL HVIDASDEEW WQARRVHSDS ETDDIGFIPS KRRVERREWS RLKAKDWGSS SGSQGREDSV LSYETVTQME VHYPARPIIL GPTKDRANDD LLSEFPDKFG SCVPHTRPK REYEIDGRDY HFVSSREKME KDIQAHKFIE AGQYNSHLYG TSVQSVREVA EQGKHCILDV SANAVRRLQA AHLHPAIFI RPRSLNVLE INKRITEEQA RKAFDRATKL EQEFTECFSA IVEGDSFEEI YHKVKRVIED

image



Background

Product Name: Recombinant Human DLG4 Protein
Catalog #: PEH4506



Alternative Names

Disks large homolog 4 (Postsynaptic density protein 95) (PSD-95) (Synapse-associated protein 90) (SAP-90) (SAP90)

Note

For Research Use Only , Not for Diagnostic Use