


Product Name: Recombinant Human HSPG2 Protein
Catalog #: PEH4508



Summary

Name	PGBM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P98160
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HGIAESDWHL EGSGGNDAPG QYGAYFHDDG FLAFPGHVFS RSLPEVPETI ELEVRTSTAS GLLLWQGVEV GEAGQGKDFI SLGLQDGHLV FRYQLGSGEA RLVSEDPIND GEWHRVTALR EGRRGSIQVD GEELVSGRSP GPNVAVNAKG SVYIGGAPDV ATLTGGRFSS GITGCVKNLV LHSARPGAPP PQPLDLQHRA

image



Background

Alternative Names Basement membrane-specific heparan sulfate proteoglycan core protein (HSPG) (Perlecan) (PLC) [Cleaved into: Endorepellin; LG3 peptide]

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human HSPG2 Protein
Catalog #: PEH4508

