



Product Name: Recombinant Human TGFBR3 Protein
Catalog #: PEH4512



Summary

Name	TGBR3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q03167
Host	E.coli
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TGVSRLFLVS EGSVVQFSSA NFSLTAETEE RNFPHGNEHL LNWARKEYGA VTSFTELKIA RNIYIKVGED QVFPPKCNIG KNFLSLNYLA EYLQPKAAEG CVMSSQPQNE EVHIIELITP NSNPYSAFQV DITIDIRPSQ EDLEVVKNLI LILKCKKSVN WVIKSFDVKG SLKIIAPNSI GFGKESERSM TMTKSIRDDI PSTQGNLVKW ALDNGYSPIT SYTMAPVANR FHLRLNNAE EMGDEEVHTI PPELRILLDP GALPALQNPP IRGGEGQNGG LPFPFPDISR RVWNEEGEDG LPRPKDPVIP SIQLFPGLRE PEEVQGSVDI

image



Background

Alternative Names	Transforming growth factor beta receptor type 3 (TGF-beta receptor type 3) (TGFR-3) (Betaglycan) (Transforming growth factor beta receptor III) (TGF-
--------------------------	---

Product Name: Recombinant Human TGFBR3 Protein
Catalog #: PEH4512



beta receptor type III)

Note

For Research Use Only , Not for Diagnostic Use