


Product Name: Recombinant Human MVK Protein
Catalog #: PEH4514



Summary

Name	KIME_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q03426
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PGKVLHGEH AVVHGKVALA VSLNLRNFLR LQPHSNGKVD LSLPNIGIKR AWDVARLQSL DTSFLEQGDV TTPTSEQVEK LKEVAGLPDD CAVTERLAVL AFLLYLSIC RKQRALPSLD IVVWSELPPG AGLGSSAAYS VCLAAALLTV CEEIPNPLKD GDCVNRWTKD DLELINKWAF QGERMIHGNP SGVDNAVSTW GGALRYHQGK ISSLRSPAL QILLTNTKVP RNTRALVAGV RNRLKFPDI VAPLLTSIDA ISLECERVLG EMGEAPAPEQ YLVLEELIDM NQHHLNALGV GHASLDQLCQ VTRARGLHSE LTGAGGGGCG ITLLKPGLEQ PEVEATKQAL TSCGFDCLST SIGAPGVSIH SATSLDSRVQ

image



Background

Alternative Names	Mevalonate kinase (MK) (EC 2.7.1.36)
--------------------------	--------------------------------------

Product Name: Recombinant Human MVK Protein
Catalog #: PEH4514



Note

For Research Use Only , Not for Diagnostic Use