


Product Name: Recombinant Human TAP1 Protein
Catalog #: PEH4515



Summary

Name	TAP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q03518
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YLDRTPRCPP SGLLTPLHLE GLVQFQDVSF AYPNRPDVLV LQGLTFTLRP GEVTALVGPN GSGKSTVAAL LQNLQPTGG QLLDGKPLP QYEHRYLHRQ VAAVGQEPQV FGRSLQENIA YGLTQKPTME EITAAAVKSG AHSFISGLPQ GYDTEVDEAG SQLSGGQRQA VALARALIRK PCVLILDDAT SALDANSQLQ VEQLLYESPE RYSRSVLLIT QHLSLVEQAD HILFLEGGAI REGGTHQQLM EKKGCYWAMV

image



Background

Alternative Names	Antigen peptide transporter 1 (APT1) (ATP-binding cassette sub-family B member 2) (Peptide supply factor 1) (Peptide transporter PSF1) (PSF-1) (Peptide transporter TAP1) (Peptide transporter involved in antigen
--------------------------	--

Product Name: Recombinant Human TAP1 Protein
Catalog #: PEH4515



processing 1) (Really interesting new gene 4 protein)

Note

For Research Use Only , Not for Diagnostic Use