


**Product Name: Recombinant Human SULT2A1 Protein**  
**Catalog #: PEH4517**



## Summary

<b>Name</b>	ST2A1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-SUMO & N-His
<b>Accession #</b>	Q06520
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	43 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	ETLRKVRDEF VIRDEDVIL TYPKSGTNWL AEILCLMHSK GDAKWIQSVP IWERSPWVES EIGYTALSET ESPRLFSSHL PIQLFPKSFF SSKAKVIYLM RNPRDVLVSG YFFWKNMKFI KKPKSWE EYF EWFCQGTVLY GSWFDHIHG W MPMREEKNFL LLSYEELKQD TGRTIEKICQ FLGKTLEPEE LNLILKNSSF QSMKENKMSN YSLLSVDYVV DKAQLLRKGV SGDWKNHFTV AQAEDFDKLF QEKMADLPRE

## image



## Background

<b>Alternative Names</b>	Sulfotransferase 2A1 (ST2A1) (EC 2.8.2.2) (Bile salt sulfotransferase) (EC 2.8.2.14) (Dehydroepiandrosterone sulfotransferase) (DHEA-ST) (DHEA-ST8) (Hydroxysteroid Sulfotransferase) (HST) (ST2) (SULT2A3)
--------------------------	---

**Product Name: Recombinant Human SULT2A1 Protein**  
**Catalog #: PEH4517**

---



**Note**

For Research Use Only , Not for Diagnostic Use