


Product Name: Recombinant Human PLA2G7 Protein
Catalog #: PEH4524



Summary

Name	PAFA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13093
Host	E.coli
Species	Human
Molecular Mass	49 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FDWQYINPVA HMKSSAWVNK IQVLMAAASF GQTKIPRGNG PYSVGCTDLM FDHTNKGTF L RLYYPSQDND RLDTLWIPNK EYFWGLSKFL GTHWLMGNIL RLLFGSMTTP ANWNSPLRPG EKYPLVVFH GLGAFRTLYS AIGIDLASHG FIVA AVEHRD RSATATYYFK DQSA AEIGDK SWLYLRTLKQ EEETHIRNEQ VRQRAKECSQ ALSLILDIDH GKPVKNALDL KFDMEQLKDS IDREKIAVIG HSFGGATVIQ TLSEDQRFRC GIALDAWMFP LGDEVYSRIP QPLFFINSEY FQYPANIIKM KKCYSPPDKER KMITIRGSVH QNFADFTFAT GKIIGHMLKL KGDIDSNVAI DLSNKASLAF LQKHLGLHKD FDQWDCLIEG DDENLIPGTN INTTNQHIML QNSSGIEKYN

image



Background

Product Name: Recombinant Human PLA2G7 Protein
Catalog #: PEH4524



Alternative Names

Platelet-activating factor acetylhydrolase (PAF acetylhydrolase) (EC 3.1.1.47)
(1-alkyl-2-acetylglycerophosphocholine esterase) (2-acetyl-1-
alkylglycerophosphocholine esterase) (Group-VIIA phospholipase A2) (gVIIA-
PLA2) (LDL-associated phospholipase A2) (LDL-PLA(2)) (PAF 2-acylhydrolase)

Note

For Research Use Only , Not for Diagnostic Use