


Product Name: Recombinant Human PCOLCE Protein
Catalog #: PEH4528



Summary

Name	PCOC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15113
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GRLEKAQGTL TTPNWPESDY PPGISCSWHI IAPPDQVIAL TFEKFDLEPD TYCRYDSVSV FNGAVSDDSR RLGKFCGDAV PGSISSEGNE LLVQFVSDLS VTADGFSASY KTLPRGTAK GQGGPKRGT EPKVKLPPKS QPPEKTEESP

image



Background

Alternative Names	Procollagen C-endopeptidase enhancer 1 (Procollagen COOH-terminal proteinase enhancer 1) (PCPE-1) (Procollagen C-proteinase enhancer 1) (Type 1 procollagen C-proteinase enhancer protein) (Type I procollagen COOH-terminal proteinase enhancer)
--------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Note

Product Name: Recombinant Human PCOLCE Protein
Catalog #: PEH4528



For Research Use Only , Not for Diagnostic Use