


Product Name: Recombinant Human PMVK Protein
Catalog #: PEH4529



Summary

| | |
|------------------------------|---|
| Name | PMVK_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q15126 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 23 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | MAPLGGAPRL VLLFSGKRKS GKDFVTEALQ SRLGADVCAV LRLSGPLKEQ YAQEHLNLFQ RLLDTSTYKE AFRKDMIRWG EEKQADPGF FCRKIVEGIS QPIWLVS DTR RVSDIQWFRE AYGAVTQTVR VVALEQSRQQ RGWVFTPGVD DAESECGLDN FGDFDWVIEN HGVEQRLEE Q LENLIEFIRS RL |

image



Background

Alternative Names Phosphomevalonate kinase (PMKase) (hPMK) (EC 2.7.4.2)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PMVK Protein
Catalog #: PEH4529

