


Product Name: Recombinant Human PON3 Protein
Catalog #: PEH4530



Summary

Name	PON3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15166
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SREVEPVEPE NCHLIEELES GSEDIDILPS GLAFISSGLK YPGMPNFAPD EPGKIFLMDL NEQNPRAQAL EISGGFDKEL FNPHGISIFI DKDNTVYLYV VNHPHMKSTV EIFKFEEQQR SLVYLKTIKH

image



Background

Alternative Names Serum paraoxonase/lactonase 3 (EC 3.1.1.2) (EC 3.1.1.81) (EC 3.1.8.1)

Note

For Research Use Only , Not for Diagnostic Use