


Product Name: Recombinant Human TOMM20 Protein
Catalog #: PEH4531



Summary

Name	TOM20_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15388
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DPNFKNRLRE RRKKQKLAKE RAGLSKLPDL KDAEAVQKFF LEEIQLGEEL LAQGEYEKGV DHLTNAIAVC GQPQQLLQVL QQTLPVPVFQ MLLTKLPTIS QRIVSAQSLA

image



Background

Alternative Names	Mitochondrial import receptor subunit TOM20 homolog (Mitochondrial 20 kDa outer membrane protein) (Outer mitochondrial membrane receptor Tom20)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TOMM20 Protein
Catalog #: PEH4531

