


Product Name: Recombinant Human FCN2 Protein
Catalog #: PEH4533



Summary

Name	FCN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q15485
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPRTCKDLLD RGHFLSGWHT IYLPDCRPLT VLCDMDTDGG GWTVFQRRVD GSVDFYRDWA TYKQGFGSRL GEFWLGNDNI HALTAQGTSE LRVDLVDFED NYQFAKYRSF KVADEAEKYN LVLGAFVEGS AGDSLTFHNN QSFSTKDQDN DLNTGNCAVM FQGAWWYKNC HVSNLNGRYL RGTHGSFANG INWKSGKGYN YSYKVSEMKV

image



Background

Alternative Names	Ficolin-2 (37 kDa elastin-binding protein) (Collagen/fibrinogen domain-containing protein 2) (EBP-37) (Ficolin-B) (Ficolin-beta) (Hucolin) (L-ficolin) (Serum lectin p35)
--------------------------	---

Product Name: Recombinant Human FCN2 Protein
Catalog #: PEH4533



Note

For Research Use Only , Not for Diagnostic Use