


Product Name: Recombinant Human IGFBP7 Protein
Catalog #: PEH4537



Summary

Name	IBP7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16270
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SIVTPPKDIW NVTGAQVYLS CEVIGIPTPV LIWNVKVRGH YGVQRTELLP GDRDNLAIQT RGGPEKHEVT GWVLVSPLSK EDAGEYECHA SNSQGQASAS AKITVVDALH EIPVKKGEGA

image



Background

Alternative Names Insulin-like growth factor-binding protein 7 (IBP-7) (IGF-binding protein 7) (IGFBP-7) (IGFBP-rP1) (MAC25 protein) (PGI2-stimulating factor) (Prostacyclin-stimulating factor) (Tumor-derived adhesion factor) (TAF)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human IGFBP7 Protein
Catalog #: PEH4537

