


Product Name: Recombinant Human KYNU protein
Catalog #: PEH4542



Summary

Name	KYNU_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q16719
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DTVQRIAAEL KCHPTDERVA LHLDEEDKLR HFRECFYIPK IQDLPPVDLS LVNKDENAIY FLGNSLGLQP KMKVKTYLEEE LDKWAKIAAY GHEVGKRPWI TGDESIVGLM KDIVGANKEKE IALMNALTVN LHLMLSFYK PTPKRYKILL EAKAFPSDHY AIESQLQLHG LNIEESMRMI KPREGEETLR IEDILEVIEK

image



Background

Alternative Names Kynureninase (EC 3.7.1.3) (L-kynurenine hydrolase)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human KYNU protein
Catalog #: PEH4542

