


Product Name: Recombinant Human GSTA3 Protein
Catalog #: PEH4543



Summary

| | |
|------------------------------|--|
| Name | GSTA3_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q16772 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 26 kDa |
| Buffer | Lyophilized from Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | MAGKPKLHYF NGRGRMEPIR WLLAAAGVEF EEKFIGSAED LGKLRNDGSL MFQQVPMVEI DGMKLVQTRA ILNYIASKYN LYGKDIKERA LIDMYTEGMA DLNEMILLLP LCRPEEKDAK IALIKEKTKS RYFPAFEKVL QSHGQDYLVG NKLSRADISL VELLYYVEEL DSSLISNFPL LKALKTRISN LPTVKKFLQP GSPRKPPADA KALEEARKIF RF |

image



Background

Alternative Names Glutathione S-transferase A3 (EC 2.5.1.18) (GST class-alpha member 3)
(Glutathione S-transferase A3-3)

Note

Product Name: Recombinant Human GSTA3 Protein
Catalog #: PEH4543



For Research Use Only , Not for Diagnostic Use