


Product Name: Recombinant Human HAGH Protein
Catalog #: PEH4546



Summary

Name	GLO2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q16775
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AALGAACARR GLGPALLGVF CHTDLRKNLT VDEGTMKVEV LPALTDNYMY LVIDDETKEA AIVDPVQPQK VVDAARKHGV KLTTVLTTHH HWDHAGGNEK LVKLESGLKV YGGDDRIGAL THKITHLSTL QVGSLNVKCL ATPCHTSGHI CYFVSKPGGS EPPAVFTGDT Lfvagcgkfy EGTADeMCKA LLEVELGRLPP DTRVYCGHEY TINNLKFARH VEPGNAAIRE KLAwAKEKYS IGEPTVPSTL AEEFTYNPFM RVREKTVQQH AGETDPVTTM RAVRREKdqf KMPRD

image



Background

Alternative Names Hydroxyacylglutathione hydrolase, mitochondrial (EC 3.1.2.6) (Glyoxalase II) (Glx II)

Product Name: Recombinant Human HAGH Protein
Catalog #: PEH4546



Note

For Research Use Only , Not for Diagnostic Use