


Summary

Name	PPM1L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5SGD2
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSRDAVKMVK GKVAEIMQND RLGGLDVLEA EFSKTWEFKN HNVAVYSIQG RRDHMEDRFE VLTDLANKTH PSIFGIFDGH GGETAAEYVK SRLPEALKQH LQDYEKDKEN SVLSYQTILE QQILSIDREM LEKLTVSUDE AGTTCLIAL SDKDLTVANV GDSRGVLC DK DGNAIPLSHD HKPYQLKERK RIKRAGGFIS FNGSWRVQGI LAMSRSLGDY PLKLNLVVIP DPDILTFDLD KLQPEFMILA SDGLWDAFSN EEA VRFIKER LDEPHFGAKS IVLQSFYRGC PDNITVMVVK

image



Background

Alternative Names	Protein phosphatase 1L (EC 3.1.3.16) (Protein phosphatase 1-like) (Protein phosphatase 2C isoform epsilon) (PP2C-epsilon)
--------------------------	---

Product Name: Recombinant Human PPM1L Protein
Catalog #: PEH4548



Note

For Research Use Only , Not for Diagnostic Use