


Product Name: Recombinant Human VWA1 Protein
Catalog #: PEH4553



Summary

Name	VWA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6PCB0
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PRGDLMFLLD SSASVSHYEF SRVREFVGQL VAPLPLGTGA LRASLVHVGS RPYTEFPFGQ HSSGEAAQDA VRASAQRMGD THTGLALVYA KEQLFAEASG ARPGVPKVLV WVTGGSSDP VGPPMQELKD LGVTVFIVST GRGNFLELSA AASAPAEKHL HFVDVDDLHI IVQELRGSIL DAMRPQQLHA TEITSSGFRL AWPPLLADS GYYVLELVPS AQPGAARRQQ LPGNATDWIW AGLDPDPTDYG VALVPESNVR LLRPQILRVR TRPGEAGPGA SGPESGAGPA PTQLAALPAP

image



Background

Alternative Names von Willebrand factor A domain-containing protein 1

Product Name: Recombinant Human VWA1 Protein
Catalog #: PEH4553



Note

For Research Use Only , Not for Diagnostic Use