


Product Name: Recombinant Human BTN2A1 Protein
Catalog #: PEH4560



Summary

Name	BT2A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7KYR7
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLQKEKKILS GEKEFERETR EIALKELEKE RVQKEEELQV KEKLQEELRW RRTFLHAVDV VLDPDTAHPD LFLSEDRRSV RRCPFRLHGE SVPDNPERFD SQPCVLGRES FASGKHYWEV EVENVIEWTV GVCRRDSVERK GEVLLIPQNG FWTLEMHKGQ YRAVSSPDRI LPLKESLCRV GVFLDYEAGD VSFYNMRDRS HIYTCPRSAF SVPVRPFRL GCEDSPIFIC PALTGANGVT VP EEGTLHR

image



Background

Alternative Names Butyrophilin subfamily 2 member A1

Note

Product Name: Recombinant Human BTN2A1 Protein
Catalog #: PEH4560



For Research Use Only , Not for Diagnostic Use