


Product Name: Recombinant Human MAEA Protein
Catalog #: PEH4561



Summary

Name	MAEA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q7L5Y9
Host	E.coli
Species	Human
Molecular Mass	71 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LKRKAVESIQ AEDESAKLCK RRIEHLKEHS SDQPAAASVW KRKRMDRMMV EHLRCGYYN TAVKLARQSG IEDLVNIEMF LTAKEVEESL ERRETATCLA WCHDNKSRLR KMKSCLEFSL RIQEFIELIR QNKRLDAVRH ARKHFSQAEG SQLDEVRQAM GMLAFPPDTH ISPYKDLLDP ARWRMLIQQF RYDNYRLHQL GNNSVFTLTL QAGLSAIKTP QCYKEDGSSK SPDCPVCSRS

image



Background

Alternative Names	E3 ubiquitin-protein transferase MAEA (EC 2.3.2.27) (Cell proliferation-inducing gene 5 protein) (Erythroblast macrophage protein) (Human lung cancer oncogene 10 protein) (HLC-10) (Macrophage erythroblast attacher) (P44EMLP)
--------------------------	--

Product Name: Recombinant Human MAEA Protein
Catalog #: PEH4561



Note

For Research Use Only , Not for Diagnostic Use