


Product Name: Recombinant Human CNTN4 Protein
Catalog #: PEH4563



Summary

Name	CNTN4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8I WV2
Host	E.coli
Species	Human
Molecular Mass	54 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GEYEPKIEVQ FPETVPTAKG ATVKLECFAL GNPVPTIWR RADGKPIARK ARRHKSNGIL EIPNFQQEDA GLYECVAENS RGKNVARGQL TFYAQPNWVQ KINDIHVAME ENVFWECKAN GRPKPTYKWL KNGEPLLTRD RIQIEQGTLN ITIVNLS DAG MYQCLAENKH GVIFSNAELS VIAVGPDFSR TLLKRVTLVK VGGEVVIECK PKASPKPVYT WKKGRDILKE NERITISEDG NLRINVTKS DAGSYTCIAT NHFGTASSTG NLVVKDPTRV MVPPSSMDVT VGESIVLPCQ VTHDHSLDIV FTWSFNGLHI DFDRDGDHFE RVGGQDSAGD LMIRNIQLKH AGKYVCMVQT SVDRLSAAAD LIVRGPPGPP EAVTIDEITD TTAQLSWRPG PDNHSPITMY VIQARTPFSV GWQAVSTVPE LIDGKTFTAT VVGLNPWVEY EFRTVAANVI GIGEPSRPSE

image



Background

Product Name: Recombinant Human CNTN4 Protein
Catalog #: PEH4563



Alternative Names

Contactin-4 (Brain-derived immunoglobulin superfamily protein 2) (BIG-2)

Note

For Research Use Only , Not for Diagnostic Use