



Product Name: Recombinant Human AGR3 Protein
Catalog #: PEH4572



Summary

| | |
|------------------------------|---|
| Name | AGR3_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | N-MBP & C-His |
| Accession # | Q8TD06 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 59 kDa |
| Buffer | Lyophilized from Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | IAIKKEKRPP QTLSRGWGDD ITWVQTYEEG LFYAQKSKKP LMVIHHLEDC QYSQALKKVF AQNEEQEMA QNKFIMLNLM HETTDKNLSP DGQYVPRIMF VDPSLTVRAD IAGRYSNRLY TYEPRDLPLL IENMKKALRL IQSEL |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Anterior gradient protein 3 (AG-3) (AG3) (hAG-3) (Anterior gradient 3 homolog) (Breast cancer membrane protein 11) (Protein disulfide isomerase family A, member 18) |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AGR3 Protein
Catalog #: PEH4572

