


Product Name: Recombinant Human AGR3 Protein
Catalog #: PEH4573



Summary

Name	AGR3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q8TD06
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IAIKKEKRPP QTLSRGWGDD ITWVQTYEEG LFYAQKSKKP LMVIHHLEDC QYSQALKKVF AQNEEIQEMA QNKFIMLNLM HETTDKNLSP DGQYVPRIMF VDPSLTVRAD IAGRYSNRLY TYEPRDLPLL IENMKKALRL IQSEL

image



Background

Alternative Names	Anterior gradient protein 3 (AG-3) (AG3) (hAG-3) (Anterior gradient 3 homolog) (Breast cancer membrane protein 11) (Protein disulfide isomerase family A, member 18)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AGR3 Protein
Catalog #: PEH4573

