

Product Name: Recombinant Human TSEN15 Protein
Catalog #: PEH4574



Summary

Name	SEN15_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WW01
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEERGDSEPT PGCSGLPGGG VRFGDGGGA PSWAPEDAWM GTHPKYLEMM ELDIGDATQV YVAFVYLDL MESKSWHEVN CVGLPELQLI CLVGTEIEGE GLQTVVPTPI TASLSHNRIR EILKASRKLQ GDPDLPMSFT LAIVESDSTI VYYKLTDFGM LPDPQNISLR R

image



Background

Alternative Names tRNA-splicing endonuclease subunit Sen15 (SEN15 homolog) (HsSEN15)
(tRNA-intron endonuclease Sen15)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TSEN15 Protein
Catalog #: PEH4574

