


Product Name: Recombinant Human RTN4IP1 Protein
Catalog #: PEH4575



Summary

Name	RT4I1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8WWV3
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YGKNEVLRFT QNMMMPIIHY PNEVIVKVHA ASVNPIDVNM RSGYGATALN MKRDPLHVKI KGEEFPLTLG RDVSGVVMCE GLDVKYFKPG DEVWAAVPPW KQGTLSSEFVV VSGNEVSHKP

image



Background

Alternative Names	Reticulon-4-interacting protein 1, mitochondrial (NOGO-interacting mitochondrial protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use