


Product Name: Recombinant Human ATG4A Protein
Catalog #: PEH4576



Summary

Name	ATG4A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WYN0
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EYPDTDELVW ILGKQHLLKT ESKLLSDIS ARLWFTYRRK FSPIGGTGPS SDAGWGCMRLR CGQMMLAQAL ICRHLGRDWS WEKQKEQPKE YQRILQCFLD RKDCCYSIHQ MAQMGVGEGK SIGEWFGPNT VAQVLKKLAL FDEWNSLAVY VSMDN TVVIE DIKKMCRVLP LSADTAGDRP PDSLTASNQS KGTSAYCSAW KPLLLIVPLR LGINQINPVY VDAFKECFKM PQSLGALGGK PNNAYYFIGF LGDELIFLDP HTTQTFVDTE

image



Background

Alternative Names	Cysteine protease ATG4A (EC 3.4.22.-) (AUT-like 2 cysteine endopeptidase) (Autophagin-2) (Autophagy-related cysteine endopeptidase 2) (Autophagy-related protein 4 homolog A) (hAPG4A)
--------------------------	--

Product Name: Recombinant Human ATG4A Protein
Catalog #: PEH4576



Note

For Research Use Only , Not for Diagnostic Use