


Product Name: Recombinant Human ST8SIA1 Protein
Catalog #: PEH4578



Summary

Name	SIA8A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q92185
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LPNEKEIVQG VLQQGTAWRR NQTAARAFRK QMEDCCDPAH LFAMTKMNSP MGKSMWYDGE FLYSFTIDNS TYSLFPQATP FQLPLKCAV VGNGGILKKS GCGRQIDEAN FVMRCNLPL SSEYTKDVGS

image



Background

Alternative Names	Alpha-N-acetylneuraminide alpha-2,8-sialyltransferase (EC 2.4.99.8) (Alpha-2,8-sialyltransferase 8A) (Ganglioside GD3 synthase) (Ganglioside GT3 synthase) (Sialyltransferase 8A) (SIAT8-A) (Sialyltransferase St8Sia I) (ST8Sial)
--------------------------	--

Note

Product Name: Recombinant Human ST8SIA1 Protein
Catalog #: PEH4578



For Research Use Only , Not for Diagnostic Use