


Product Name: Recombinant Human FAM3C Protein
Catalog #: PEH4580



Summary

Name	FAM3C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92520
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QVFEIKMDAS LGNLFARSAL DTAARSTKPP RYKCGISKAC PEKHFAFKMA SGAANVVGPK ICLEDNVLMS GVKNNVGRGI NVALANGKTG EVLDTKYFDM WGGDVAPFIE FLKAIQDGTI VLMGTYDDGA TKLNDEARRL IADLGSTSIT NLGFRDNWVF CGGKGIKTKS PFEQHIKNNK DTNKYEGWPE VVEMEGCIPQ KQD

image



Background

Alternative Names Protein FAM3C (Interleukin-like EMT inducer)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FAM3C Protein
Catalog #: PEH4580

