


Product Name: Recombinant Human NCLN Protein
Catalog #: PEH4584



Summary

Name	NCLN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q969V3
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SRRCVLMRLL DFSYEQYQKA LRQSAGAVVI ILPRAMAAVP QDVVRQFMEI EPEMLAMETA VPVYFAVEDE ALLSIYKQTQ AASASQGSAS AAEVLLRTAT ANGFQMVTSG VQSKAVSDWL IASVEGRLTG LGGEDLPTIV IVAHYDAFGV APWLSLGADS NGSGVSVLLE LARLFSRLYT YKRTHAAYNL LFFASGGGKF NYQGTKRWLE DNL DHTDSSL LQDNVAFVLC LDTVGRGSSL HLHVSKPPRE GTLQHAFLRE LETVAHQFP EVRFSMVHKKR INLAEDVLAW

image



Background

Alternative Names Nicalin (Nicastrin-like protein)

Product Name: Recombinant Human NCLN Protein
Catalog #: PEH4584



Note

For Research Use Only , Not for Diagnostic Use