


Product Name: Recombinant Human RAB37 Protein
Catalog #: PEH4585



Summary

Name	RAB37_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96AX2
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MTGTPGAVAT RDGEAPERSP PCSPSYDLTG KVMLLGDTGV GKTCFLIQFK DGAFLSGTFI ATVGIDFRNK VVTVDGVRVK LQIWDTAGQE RFRSVTHAYY RDAQALLLY DITNKSSFDN IRAWLTEIHE YAQRDVVIML LGNKADMSSE RVIRSEDGET LAREYGV PFL ETS AKTGMNV ELAFLAIAKE LKYRAGHQAD EPSFQIRDYV ESQKKRSSCC

image



Background

Alternative Names Ras-related protein Rab-37

Note

Product Name: Recombinant Human RAB37 Protein
Catalog #: PEH4585



For Research Use Only , Not for Diagnostic Use