


Product Name: Recombinant Human RGMA Protein
Catalog #: PEH4586



Summary

Name	RGMA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q96B86
Host	E.coli
Species	Human
Molecular Mass	54 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CKILKCNSEF WSATSGSHAP ASDDTPEFCA ALRSYALCTR RTARTCRGDL AYHSAVHGIE DLMSQHNCSK DGPTSQPRLR TLPPAGDSQE RSDSPEICHY EKSFHKHSAT PNYTHCGLFG DPHLRTFTDR FQTCKVQGAW PLIDNNYLVN QVTNTPVLPG SAATATSKLT IIFKNFQECV DQKVYQAEMD ELPAAFVDGS KNGGDKHGAN SLKITEKVSG QHVEIQAKYI GTTIVVRQVG RYLTFAVRMP EEVVNAVEDW DSQGLYLCLR GCPLNQQIDF QAFHTNAEGT GARRLAAASP APTAPETFPY ETAVAKCKEK LPVEDLYYQA CVFDLLTTGD VNFTLAAYYA LEDVKMLHSN KDKLHLYERT RDLPG

image



Background

Alternative Names	Repulsive guidance molecule A (RGM domain family member A)
--------------------------	--

Product Name: Recombinant Human RGMA Protein
Catalog #: PEH4586



Note

For Research Use Only , Not for Diagnostic Use