


**Product Name: Recombinant Human RNF31 Protein**  
**Catalog #: PEH4588**



## Summary

<b>Name</b>	RNF31_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q96EP0
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	13 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	VIEQKEVPNG LRDEACGKET PAGYAGLCQA HYKEYLVSLI NAHSLDPATL YEVEELETAT ERYLHVRPQP LAGEDPPAYQ ARLLQKLTEE VPLGQSIPRR

## image



## Background

<b>Alternative Names</b>	E3 ubiquitin-protein ligase RNF31 (EC 2.3.2.31) (HOIL-1-interacting protein) (HOIP) (RING finger protein 31) (RING-type E3 ubiquitin transferase RNF31) (Zinc in-between-RING-finger ubiquitin-associated domain protein)
--------------------------	---

## Note

For Research Use Only , Not for Diagnostic Use