


Product Name: Recombinant Human SAT2 Protein
Catalog #: PEH4589



Summary

Name	SAT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96F10
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MASVRIREAK EGDCGDILRL IRELAEFEKL SDQVKISEEA LRADGFGDNP FYHCLVAEIL PAPGKLLGPC VVGYGIIYFI YSTWKGRTIY LEDIYVMPEY RGQGIGSKII KKVAEVALDK GCSQFRLAVL DWNQRAMDLY KALGAQDLTE AEGWHFFCFQ GEATRKLAKG

image



Background

Alternative Names Thialysine N-epsilon-acetyltransferase (EC 2.3.1.-) (Diamine acetyltransferase 2) (EC 2.3.1.57) (Spermidine/spermine N(1)-acetyltransferase 2) (SSAT-2)

Note

Product Name: Recombinant Human SAT2 Protein
Catalog #: PEH4589



For Research Use Only , Not for Diagnostic Use