


Product Name: Recombinant Human ERMAP Protein
Catalog #: PEH4593



Summary

Name	ERMAP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96PL5
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VTEVDNLLSD HAKEKGKLHK AVKKLRSELK LKRAAANSWG RRARLHFVAV TLDPDTAHPK LILSEDQRCV RLGDRRQVPV DNPQRFDVAV SILGSEYFTT GCHYWEVYVG DKTKWILGVC SESVSRKGKV TASPANGHWL LRQSRGNEYE ALTSPQTSFR LKEPPRCVGI FLDYEAGVIS FYNVTNKSHI FTFTHNFSGP LRPFEPCLH DGGKNTAPLV ICSELHKSEE SIVPRPEGKG HANGDVSLKV NSSLLPPKAP ELKDIILSLP

image



Background

Alternative Names	Erythroid membrane-associated protein (hERMAP) (Radin blood group antigen) (Scianna blood group antigen)
--------------------------	--

Product Name: Recombinant Human ERMAP Protein
Catalog #: PEH4593



Note

For Research Use Only , Not for Diagnostic Use