


Product Name: Recombinant Human SORCS2 Protein
Catalog #: PEH4594



Summary

Name	SORC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q96PQ0
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGGVDMQQSQ VQLQCPLTPP RGLQVSIQGE AVAVRPGEDV LFVVRQEQGD VLTTYQVDL GDGFKAMYVN LTLTGPIRH RYESPGIYRV SVRAENTAGH DEAVLFVQVN SPLQALYLEV VPVIGLNQEV NLTAVLLPLN PNLTVFYWWI GHSLQPLLSL DNSVTTRFSD TGDVRVTVQA ACGNSVLQDS RVLRLDQFQ VMPLQFSKEL DAYNPNTPEW REDVGLVVTR LLSKETSV PQ ELLVTVVKPG LPTLADLYVL LPPRPTRKR SLSSDKRLAA IQQVLNAQKI

image



Background

Alternative Names	VPS10 domain-containing receptor SorCS2 [Cleaved into: SorCS2 122 kDa chain; SorCS2 104 kDa chain; SorCS2 18 kDa chain]
--------------------------	---

Product Name: Recombinant Human SORCS2 Protein
Catalog #: PEH4594



Note

For Research Use Only , Not for Diagnostic Use