


Product Name: Recombinant Human COQ7 Protein
Catalog #: PEH4596



Summary

Name	COQ7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q99807
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDHLKKFNEL MVTRVRPTV LMPLWNVLGF ALGAGTALLG KEGAMACTVA VEESIAHHYN NQIRTLMEED PEKYEELLQL IKKFRDEELE HHDIGLDHDA

image



Background

Alternative Names	5-demethoxyubiquinone hydroxylase, mitochondrial (DMQ hydroxylase) (EC 1.14.99.60) (Timing protein clk-1 homolog) (Ubiquinone biosynthesis monooxygenase COQ7)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use