


Product Name: Recombinant Human SH3GL3 Protein
Catalog #: PEH4597



Summary

Name	SH3G3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99963
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSVAGLKKQF HKASQLFSEK ISGAEGTKLD DEFLDMERKI DVTNKVVAEI LSKTTEYLQP NPAYRAKLGM LNTVSKIRGQ VKTTGYPQTE GLLGDCMLKY GKELGEDSTF GNALIEVGES MKLMAEVKDS LDINVKQTFI DPLQLLQDKD LKEIGHHLKK LEGRRLDYDY KKKRVGKIPD EEVRQAVEKF EESKELAERS MFNFLENDVE QVSQLAVFIE AALDYHRQST EILQELQSKL QMRISAASSV PRREYKPRPV KRSSSELNGV STTSVVKTTG SNIPMDQPCC RGLYDFEPEN QGELGFKEGD IITLTNQIDE NWYEGMIHGE SGFFPINYVE VIVPLPQ

image



Background

Alternative Names	Endophilin-A3 (EEN-B2) (Endophilin-3) (SH3 domain protein 2C) (SH3 domain-containing GRB2-like protein 3)
--------------------------	---

Product Name: Recombinant Human SH3GL3 Protein
Catalog #: PEH4597



Note

For Research Use Only , Not for Diagnostic Use