


Product Name: Recombinant Human HSD17B14 Protein
Catalog #: PEH4598



Summary

Name	DHB14_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BPX1
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VVVVTGGGRG IGAGIVRAFV NSGARVVICD KDESGGRALE QELPGAVFIL CDVTQEDDVK TLVSETIRRF GRLCDVVNNA GHHPPPQRPE ETSAQGFRQL LELNLLGTYT LTKLALPYLR KSQGNVINIS SLVGAIGQAQ AVPYVATKGA VTAMTKALAL DESPYGVRVN CISP GNIWTP LWEEAALMP DPRATIREGM LAQPLGRMGQ PAEVGAAAVF LASEANFCTG IELLVTGGAE LGYGCKASRS

image



Background

Alternative Names	17-beta-hydroxysteroid dehydrogenase 14 (17-beta-HSD 14) (EC 1.1.1.62) (17-beta-hydroxysteroid dehydrogenase DHRS10) (Dehydrogenase/reductase SDR family member 10) (Retinal short-chain dehydrogenase/reductase retSDR3) (Short chain dehydrogenase/reductase
--------------------------	---



Product Name: Recombinant Human HSD17B14 Protein
Catalog #: PEH4598

family 47C member 1)

Note

For Research Use Only , Not for Diagnostic Use