


Product Name: Recombinant Human THTPA Protein
Catalog #: PEH4601



Summary

Name	THTPA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BU02
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAQGLIEVER KFLPGPGTEE RLQELGGTLE YRVTFRDITY DPELSLMQA DHWLRRREDS GWELKCPGAA GVLGPHTEYK ELTAEPTIVA QLCKVLRADG LGAGDVA AVL GPLGLQEVAS FVTKRSAWKL VLLGADEEEP QLRVDLDTAD FGYAVGEVEA LVHEEA EVPT ALEKIHRLSS MLGVPAQETA PAKLIVYLQR FRPDYQRLLEVNSSRERPQ ETEDPDHCLG

image



Background

Alternative Names Thiamine-triphosphatase (ThTPase) (EC 3.6.1.28)

Note

Product Name: Recombinant Human THTPA Protein
Catalog #: PEH4601



For Research Use Only , Not for Diagnostic Use